

compute function:

dt-record:

1. **sml**=quantity of suppressor when addition completed

2. **sv**=total volume of suppressor added when addition completed

3. **sc**=**sml**/**sv**

4. **vv**=base volume

5. compute x value of output curve at y=**evaluation ratio** using cubic spline interpolation. let result be **sv1**.

6. calibration factor **z**=**sc**/(1+**vv**/**sv1**);

$$Z = \frac{C_S}{1 + \frac{V_b}{V_S(R)}}$$

dt-analysis:

1.vv=base volume

2. compute x value of output curve at y=evaluation ratio using cubic spline interpolation. let result be spv.

3.if calibration factor z is given, sample suppressor concentration  $spc=z*(1+vv/spv)$ .

4.if calibration factor z not given, compute z from another dt-record file with evaluation ratio in this setting. sample suppressor concentration  $spc=z*(1+vv/spv)$ .

$$C_{sp} = z \left( 1 + \frac{V_b}{V_{sp}(R)} \right)$$

lat-record:

1.compute the cubic spline interpolation of output curve.

2.compute the expression of derivative of cubic spline interpolation.

3.compute the last intersection of derivative curve and line y=thres. default value of thres is -0.05. the x value of result is intercept concentration.

lat-analysis:

1.spv=total volume of solution when sample addition completed.

2.spv0=total volume of sample added.

3.intercept Q itq=charge value when suppressor addition completed

4.select data points of accelerator addition from output curve and do linear fitting. result is  $y=k*x+b$ .

5.sample accelerator concentration  $spc=-(itq-b)/k*spv/spv0$ .

$$C_{sp} = \frac{V_{total}(b - Q_{it})}{V_{sp}k}$$

rc-record:

1.Q=evaluation ratio\*charge of last suppressor-accelerator addition step.

2.compute x value of output curve at  $y=Q$ , using cubic spline interpolation. result is leveler concentration.

rc-analysis:

1.spv=total volume of solution when sample addition completed

2.  $spv0$ =total volume of sample added
3.  $q$ =charge of last sample addition step.
4. read another rc-record file,  $q0$ =charge of last step of this rc-record file.
5. compute leveler concentration of this rc-record file with  $evaluation\ ratio=q/q0$ . result is  $lc$ .
6. sample leveler concentration  $spc=lc*spv/spv0$ .

$$C_{sp} = \frac{C_L V_{total}}{V_{sp}}$$

sar-record:

1. compute x value of the output curve at  $y=evaluation\ ratio$ , result is  $sconc$ .
2.  $ac$ =accelerator concentration of solution when suppressor-accelerator addition completed.
3.  $sc$ =suppressor concentration of solution when suppressor-accelerator addition completed.
4.  $aconc=sconc*ac/sc$ .
5. for each output curve, repeat step 1-4, and result are the data points ( $aconc,sconc$ )

6.do linear fitting to data points ( $a_{conc}$ ,  $s_{conc}$ ). result is the sa-relation  $s=k*a+b$ .

sar-analysis:

1. $vv$ =base volume

2. $sp_{vend}$ =volume of sample added

3.compute sa-relation line from another sar-record file, result line  $s=k_0*a+b_0$ .

4.select the output curve which is zero accelerator concentration from this sar-record file. let this output curve be  $snq_0$ .

5.do linear fitting to second output curve (correspond to accelerator addition). result line is  $y=k*x+b$ .

6.compute x value of first output curve (correspond to sample addition) at  $y$ =evaluation ratio using cubic spline interpolation. result is  $spv$ .

7.let  $ac=0$ .

8. $sc=k_0*ac+b_0*(vv/spv+1)$

9.compute y value of curve  $snq_0$  at  $x=sc/(vv/spv+1)$  using cubic spline interpolation. result is  $tmp$ .

10. $ac=(b-tmp)/k*(vv/spv+1)$ ;

11.repeat step 8-10  $n$  times.  $n=3$  by default.

12.sample suppressor concentration is  $sc$ , sample accelerator concentration is  $ac$ .

$$C'_{Ssp} = k_{SAR} C'_{A sp} + \frac{b_{SAR}}{G} \left( \frac{V_b}{V_{sp}(R)} + 1 \right)$$

$$C'_{A sp} = \frac{b_A - q(C_S = C'_{Ssp}, C_A = 0)}{k_A}$$

$$G = \frac{V_b}{V_{sp}} + 1$$

$$C_{Ssp} = G C'_{Ssp}$$

$$C_{Asp} = G C'_{Asp}$$

pal-record:

1. compute x value of the output curve at y=**evaluation ratio**, result is leveler concentration.

pal-analysis:

- 1.**spv**=total volume of solution when sample addition completed
- 2.**spv0**=total volume of sample added.
- 3.**lml**=quantity of leveler when leveler addition completed.
- 4.**nq**=normalized charge when sample addition completed.
5. compute leveler concentration from another pal-record file with **evaluation ratio=nq**. result is **lc**.
- 6.sample leveler concentration **spc**=(**lc\*spv-lml**)/**spv0**.

$$C_{sp} = \frac{C_L V_{total} - m_L}{V_{sp}}$$

Irt-record:

1.do linear fitting to output curve, result line is  $y=k*x+b$ , this is the regression line.

Irt-analysis:

1. $spv0$ =total volume of sample added.

2. $lml$ =quantity of leveler when leveler addition completed

3. $nq$ =normalized charge when leveler addition completed

4. $spv$ =total volume of solution when leveler addition completed

5.read another Irt-record file and compute regression line  $y=k0*x+b0$ .

6.ignore first  $n$  points of output curve and do linear fitting, result fitting line is  $y=k*x+b$ .  $n=3$  by default.

7.sample leveler concentration  $spc=((nq-b0)/k*spv-lml)/spv0$ .



$$C_{sp} = \frac{V_{total}(q - b_{std})}{kV_{sp}} - \frac{m_L}{V_{sp}}$$